



| **Contacolonie Astor 20**

Astor 20 Colony Counter

The innovative counter based on pressure impulses



The new Astor 20 Colony Counter is very easy to use. Thanks to its Petri dish holder, back-lit and sensitive to pressure, the enumeration of the microbial grown colonies on the agar surface can be performed by a simple and common marker. The Petri dish holder is back-lit by LED's: this new technology grants a better light diffusion, a greater reliability and a minimum power consumption (2 W, as a total).

Astor 20 does not require any particular marker, connected by a cable.

Technical features:

- Easy-to-read display with automatic counter and Reset key
- Acoustic signal to confirm every colony count
- Transparent and back-lit Petri dish holder, sensitive to pressure
- Every kind of commercial marker can be used: no need of a special one

- Petri dishes with diameters up to 110 mm can be used
- 90-mm Petri dish adapter included
- Minimum power consumption: 2 W
- It can be used either with 110 or 220 V

Optional accessory:

- Code no. 66131 - Benchtop lens with lamp, adjustable arm

